


### 3.4. Target Design for Advanced Alignment on an MA/BA6 Gen3

The use of Assisted and AutoAL alignment methods requires certain conditions to be met to ensure that alignment of the targets will perform correctly on the aligner. Some basic essentials for successful alignments are presented next:

#### 3.4.1. The Clear Field (TSA Auto and Assisted Alignment)



### The Field of View (FOV)

Objective	Olympus 5X	Olympus 10X	Olympus 20X
FOV	1300 x 900	650 x 460	250 x 190
Magnification	210	420	1050
Recommended Clear Area	1400 x 1000	750 x 560	350 x 290

All values in microns ( $\mu\text{m}$ )

The size of a clear area depends on two more variables:

- The **field of view**, which depends on the **objective magnification** (see size information shown in the table left) and
- The **wafer misplacement errors** of the operator, if the clear field is smaller than the field of view.

### 3.4.2. Alignment Targets for Assisted and Automatic Alignment

	<p>The key for a successful assisted and automatic alignment is the design of a <b>unique</b> target within the field of view. As a rule of thumb you can say that mask and the wafer targets should have less than 25% of the same edges.</p> <p>If the targets are not unique then the system may locate an image that is not correct and misalignments may occur.</p> <ul style="list-style-type: none"> <li><b>Target 1 in FOV <math>\neq</math> Target 2 in FOV</b></li> </ul>
--	---

### Target Size Ratio

<p>Mask (left) &amp; wafer (right) target</p> <p>Mask target aligned to wafer target</p>	<p>The targets at the left are the most commonly used alignment marks.</p> <p><b>In manual mask aligners it is essential that the mask target does not cover the wafer target completely!</b> In most cases this delta should be greater than five microns edge to edge.</p> <p>The final alignment accuracy depends greatly on the quality of mask and wafer alignment targets due to process related parameters (contrast, topography, form fidelity, etc.).</p> <ul style="list-style-type: none"> <li><math>c \geq a + 10\mu\text{m}</math> (2 x <math>5\mu\text{m}</math>)</li> <li><math>d \geq b + 10\mu\text{m}</math> (2x <math>5\mu\text{m}</math>)</li> </ul>
--	--


### Comparison Target Size to Field of View

<p><b>Mask:</b></p> <table border="1"> <thead> <tr> <th>Objective</th> <th>a [<math>\mu\text{m}</math>]</th> <th>b [<math>\mu\text{m}</math>]</th> </tr> </thead> <tbody> <tr> <td>20X TSA Only</td> <td>70</td> <td>20</td> </tr> <tr> <td>10X TSA &amp; BSA</td> <td>124</td> <td>35</td> </tr> <tr> <td>5XTSA &amp; BSA</td> <td>250</td> <td>70</td> </tr> </tbody> </table> <p>Dimensions of mask alignment targets</p> <p><b>Wafer:</b></p> <table border="1"> <thead> <tr> <th>Objective</th> <th>c [<math>\mu\text{m}</math>]</th> <th>d [<math>\mu\text{m}</math>]</th> </tr> </thead> <tbody> <tr> <td>20X TSA Only</td> <td>80</td> <td>30</td> </tr> <tr> <td>10X TSA &amp; BSA</td> <td>132</td> <td>43</td> </tr> <tr> <td>5X TSA &amp; BSA</td> <td>265</td> <td>85</td> </tr> </tbody> </table> <p>Dimensions of wafer alignment targets</p>			Objective	a [ $\mu\text{m}$ ]	b [ $\mu\text{m}$ ]	20X TSA Only	70	20	10X TSA & BSA	124	35	5XTSA & BSA	250	70	Objective	c [ $\mu\text{m}$ ]	d [ $\mu\text{m}$ ]	20X TSA Only	80	30	10X TSA & BSA	132	43	5X TSA & BSA	265	85	<p>The tables on the left show typical dimensions for mask / wafer targets based upon the objective that is used.</p>
Objective	a [ $\mu\text{m}$ ]	b [ $\mu\text{m}$ ]																									
20X TSA Only	70	20																									
10X TSA & BSA	124	35																									
5XTSA & BSA	250	70																									
Objective	c [ $\mu\text{m}$ ]	d [ $\mu\text{m}$ ]																									
20X TSA Only	80	30																									
10X TSA & BSA	132	43																									
5X TSA & BSA	265	85																									


#### 3.4.3. Working with Auto Origin

	<p><b>Auto Origin = Automatic setting of the target center</b>, which is done by the MA/BA6 Gen3 automatically after the target has been defined. When the target has been found the image is switched around by 180° and the center is calculated. This is usually more precise than a manual setting of the target center.</p> <p><b>Note!</b> For the Auto Origin function the targets need to be symmetrically.</p>
--	---

### 3.4.4. Working with Manual Origin

	<p>If the targets are <b>not symmetrically</b>, it is possible to set the origin <b>manually</b>.</p> <p>For the <b>visual alignment</b> check there are two options:</p> <ul style="list-style-type: none"> <li>• If there is enough space <b>verniers</b> (see image on left side) can get created on the mask and wafer, which can be read on the screen before exposure. This is the most precise version.</li> <li>• It is also possible to design the wafer and mask targets as described in 3.4.2. A skilled operator will be able to work with it.</li> </ul>
---	---

### 3.4.5. Proposed Target Design for Assisted and Automatic Alignment (MA/BA6 Gen3)

 <p>The diagram illustrates the target design for assisted and automatic alignment. It shows two fields of view: Left Field of View and Right Field of View. The vertical dimension is 512 μm and the horizontal dimension is 384 μm. The target design includes a Mask Target (blue square) and a Wafer Target (red square). The target sizes are 50 ... 70 μm and 5 ... 10 μm.</p>	<p><b>Proposal 1</b></p> <ul style="list-style-type: none"> <li>• <b>Independent mask / wafer targets</b></li> <li>• <b>Partial occlusion possible</b></li> </ul> <p>In the design example shown left the wafer and mask targets can be seen independently. It may occur that the wafer target is partially covered by the mask target. In this case an automatic or manual movement of the mask target may be required.</p>
---	--

<p>Left Field of View      Right Field of View</p> <p>512 μm</p> <p>384 μm</p> <p>Mask Target</p> <p>Wafer Target</p> <p>10 μm</p> <p>100 μm</p> <p>130 ... 150 μm</p>	<p><b>Proposal 2</b></p> <ul style="list-style-type: none"> <li>• Same capabilities as Proposal 1</li> <li>• Application of synthetic targets possible</li> <li>• Partial occlusion may occur</li> </ul>
<p>A = 300x400 μm, 5 μm line width</p> <p>Left Field of View      Right Field of View</p> <p>512 μm</p> <p>384 μm</p> <p>Mask Target</p> <p>Wafer Target</p> <p>B = 70 ... 100 μm</p> <p>2 ... 3 μm</p> <p>3 x 3 μm</p>	<p><b>Proposal 3</b></p> <ul style="list-style-type: none"> <li>• Proposal 3 target solution</li> </ul> <p>In this target combination no occlusion will occur. The large rectangle is used for alignment. The small squares on the mask do not cover the wafer targets (cross shaped). They have been designed for visual control.</p>

Based on Proposal 3 the following equation describes the ideal MA/BA6 Gen3 target:

$$0.5 \cdot (A-B) > \text{sqr} [(3 \cdot \sigma_{\text{PAL}})^2 + (3 \cdot \sigma_{\text{1st Exposure}})^2 + (3 \cdot \sigma_{\text{Chuck Position}})^2]$$